Safety Data Sheet

Revision: 1.06 Modified: 2020-10-22

biochrom

1: Identification of substance / mixture

1. Product Identifier Mixture

Product Name Polyamine Buffer

Product Code 80-2037-74

CAS Number

Other Names 5056336800482

IUPAC

MFCD Number EC/EINECS

REACH Number Index-No

2. Relevant identified uses of the substance or mixture and uses advised against

Laboratory use only

3. Details of the supplier of the safety data sheet

Biochrom Ltd

Unit 7, Enterprise Zone,3970 Cambridge Research Park

Waterbeach

 Cambridge
 Telephone:
 +44 (0)1223 423723

 CB25 9PE
 Fax:
 +44 (0)1223 420164

 United Kingdom
 Email:
 enquiries@biochrom.co.uk

4. Emergency telephone number

+44 (0)1223 427890 -

2. Hazards Identification

1. Classification of the substance or mixture

Non Hazardous

2. Label elements

Non Hazardous

Hazard Statements

Non Hazardous

Precautionary Phrases

Non Hazardous

3. Other Hazards

Non Hazardous

3. Composition / Information on Ingredients

2. Mixtures

Product Name		Hazards	Concentration		
Water					
	CAS Number: 7732-18-5 EC/EINECS: 231-791-2		70.6%		
Propan-2-ol					
	CAS Number: 67-63-0 EC/EINECS: 200-661-7		6.5%		
Sodium chloride					
	CAS Number: 7647-14-5 EC/EINECS: 231-598-3		3.2%		
Tri-potassium citrate					

Safety Data Sheet

Revision: 1.06 Modified: 2020-10-22

	1010dilica. 2020 10 22
CAS Number: 6100-05-6 EC/EINECS: 231-905-0	17.8%
Hydrochloric acid	
CAS Number: 7647-01-0 EC/EINECS: 231-595-7	1.8%
Phenol	
CAS Number: 108-95-2 EC/EINECS: 203-632-7	0.1%

4. First Aid Measures

1. Description of first aid measures

Skin Contact Wash immediately with plenty of soap and water.

Eye Contact Rinse opened eye for several minutes under running water.

Ingestion Do not induce vomiting.

Wash out mouth with water.

Consult a doctor.

Inhalation Supply fresh air; consult a doctor in case of complaints.

2. Most important symptoms and effects

There may be irritation and redness.

3. Indication of any immediate medical attention

No additional measures required

5. Firefighting measures

1. Extinguishing Media

Suitable Suitable extinguishing media for surrounding fire should be used.

Unsuitable N/A

2. Special Hazards arising from the substance or mixture

N/A

3. Advice for Fire Fighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear

6. Accidental Release Measures

1. Personal Precautions

Refer to section 8 of SDS for personal protection details.

2. Environmental Precautions

Do not discharge into drains or rivers.

Do not allow to enter sewers / surface or ground water.

3. Methods & Materials

Mix with sand or vermiculite.

Transfer to a closable, labelled salvage container for disposal by an appropriate method.

4. Preventing the occurrence of secondary hazards.

None

Safety Data Sheet

Revision: 1.06 Modified: 2020-10-22

7. Handling and Storage

1. Personal Precautions

Safe Handling Ensure there is sufficient ventilation of the area.

Protection against explosions and fires

No special requirements

2. Conditions for safe storage, including any incompatibilities

Managing Storage Risks Store in cool, well ventilated area.

Keep container tightly closed. Keep away from direct sunlight.

Storage Controls No special requirements

Maintaining Integrity No special requirements

Other advice No further information available.

3. Specific End Uses

The substance is supplied for Research and Development by professionals only. To be used on Biochrom Amino Acid Analyzers.

8. Exposure Controls/Personal Protection

1. Control Parameters

No Data Available

2. Exposure Controls

General protective and hygiene measures

The standard precautionary measures should be adhered to when handling

Engineering measures

Ensure there is sufficient ventilation of the area.

Eye / Face Protection

Safety Glasses.

Hand protection

Nitrile gloves.

Respiratory protection

Respiratory protection not required.

Skin protection

Protective clothing.

Other personal protection advice

No data.

9. Physical and Chemical Properties

1. Physical and Chemical Properties

Appearance	Clear liquid			
Odour	Slightly phenolic			
Odour threshold	No Data Available			
PH	5.60			
Melting point / Freezing point	No Data Available			
Initial boiling point and boiling range	No Data Available			
Flash point	No Data Available			
Evaporation rate	No Data Available			
Flammability(solid,gas)	No Data Available			
Upper/lower flammability or explosive limits	No Data Available			
Vapour pressure	No Data Available			
Vapour density	No Data Available			
Relative density	No Data Available			
Solubility(ies):	Water soluble			
Partition coefficient: n-octanol/water	No Data Available			

Safety Data Sheet

Revision: 1.06 Modified: 2020-10-22

Aut	to-ignition temperature	No Data Available
Dec	composition temperature	No Data Available
Vis	cosity	No Data Available
Exp	plosive properties	No Data Available
Oxi	idising properties	No Data Available

2. Other Information

No additional information available

10. Stability and Reactivity

1. Reactivity

No unusual reactivity

2. Stability

Stable at room temperature.

3. Possibility of Hazardous Reactions

No hazardous reactions known.

4.Conditions to Avoid

Direct Sunlight.

5. Incompatible Materials

Strong oxidizing agents.

6. Hazardous Decomposition Products

In combustion emits toxic fumes.

11. Toxicology information

1. Information

Acute Toxicity No information available

Skin corrosion/irritation No information available

Serious eye

Damage/irritation

Respiratory or skin

sensitisation

No information available

No information available

Germ Cell mutagenicity

No information available

Carcinogenicity

No information available

Reproductive toxicity

No information available

STOT-single exposure

No information available

STOT-repeated exposure

No information available

Aspiration hazard

No information available

2. Additional

No additional information available

12. Ecological Information

1. Toxicity

No information available

2. Persistence and degradability

Safety Data Sheet

Revision: 1.06 Modified: 2020-10-22

No information available

3. Bio-Accumulative Potential

No information available

4. Mobility and Soil

No information available

5. Results of PBT & vPvB assessment

No information available

6. Other adverse effects

No information available

13. Disposal Considerations

1. Waste Treatment Methods

Disposal Operations Consult state, local or national regulations for proper disposal.

Disposal of Packaging Disposal must be made according to official regulations.

14. Transport Information

Air (ICAO)

Not classified as hazardous for transport

Road (ADR)

Not classified as hazardous for transport

Sea (IMDG)

Not classified as hazardous for transport

15. Safety, health, environmental and national regulations

1. Safety, health, environmental and national regulations:

Product is not subject to any additional regulations or provisions.

2. Safety Assessment

No Chemical Safety Assessment

16. Other Information

1. Other Information:

This safety data sheet complies to the requirements of Regulation (EC) No. 1272/2008

3. Disclaimer

Safety Data Sheet

Language : English	Salety	Dala	SHEEL	N	Modified: 2020-10-22
The information provided on this SDS is correct The information given is designed only as a guand is not to be considered as a warranty or quany not be valid for such material used in com	ide for safe har Iality specificati	idling, use, on. The info	processing, storage, to the processing and the processing to the processing to the processing to the processing to the processing and the processing to the	ransportation, disportation the specific mate	osal and release rial designated and

Revision: 1.06