

# CERTIFICATE OF ANALYSIS

For

## Sodium Oxidized Buffer 1

Part Number: 80-2037-80

Batch (Control) Number: R27818

Quantity: 2 Litres

Expiry Date: 09 September 2025

The above product gave the following results on analysis:

pH = 3.346 @ 20.0°C

pH Specification = 3.33 - 3.37 @ 20.0°C

Conductivity = 15.46 mScm<sup>-1</sup> @ 20.0°C

Conductivity Specification = 13.6 - 19.2 mScm<sup>-1</sup> @ 20.0°C

Chemical QC Signature:

Date of Issue: 15 September 2022

Neil Booth

---

|   |   |                    |                        |
|---|---|--------------------|------------------------|
| Postal Address  | Telephone   | Telefax            | Email:                 |
| Biochrom Ltd<br>Unit 7, Enterprise Zone<br>3970 Cambridge Research Park<br>Beach Drive, Waterbeach<br>Cambridge<br>CB25 9PE, England. | +44 (0)1223 423723  | +44 (0)1223 420164 | support@biochrom.co.uk |
|   | Web Site: <a href="http://www.biochrom.co.uk">http://www.biochrom.co.uk</a> |                    |                        |

Registered in England No: 3526954 Registered Office: Unit 7, Enterprise Zone, 3970 Cambridge Research Park, Beach Drive, Waterbeach, Cambridge, CB25 9PE England.